



## Economy anaesthetic face masks



**Anaesthesia** ▪ Anaesthetic Face Masks

## Economy anaesthetic face masks

A range of seven colour-coded economy PVC anaesthetic face masks with a soft cushion and clear housing. Seven sizes from neonatal to extra large adult.



**Clear flexible housing and cushion**  
For good patient visibility

**Soft cushion**  
For patient comfort and air tight fit

**Reduce waste**  
Available without a hook ring

**Seven sizes available**  
From neonatal to extra large adult

Code	Code without hook ring	Description	Colour	Box Qty.
1517000	1517111	Economy, anaesthetic face mask, size 6, extra large adult, 22F	Red	25
1516000	1516111	Economy, anaesthetic face mask, size 5, large adult, 22F	Orange	30
1515000	1515111	Economy, anaesthetic face mask, size 4, medium adult, 22F	Green	40
1514000	1514111	Economy, anaesthetic face mask, size 3, small adult, 22F	Yellow	50
1513000	1513111	Economy, anaesthetic face mask, size 2, paediatric, 22F	White	35
1512000	1512111	Economy, anaesthetic face mask, size 1, infant, 15M	Black	40
1511000	1511111	Economy, anaesthetic face mask, size 0, neonate, 15M	Blue	50

Make an enquiry

IS5.6 • Issue 11 01.25



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK  
T: +44 (0)118 965 6300 info@intersurgical.com www.intersurgical.com



The manufacturer Intersurgical Ltd is certified to ISO 14001:2015, ISO 9001:2015, ISO 13485:2016 and MDSAP

**Please think before you print**  
Save energy and paper.  
If you must print this information sheet please print it double sided.